



High Performance Foams Division
 245 Woodstock Road
 Woodstock, CT 06281-1815
 Tel: 860.928.3622 / Fax: 860.928.7843
 www.rogerscorporation.com

R/bak® Cushion Mounting Materials

R/bak® Materials Proper Store and Shelf Life Guidelines

As with any materials ordered in large quantities, there is occasionally a time lapse between the purchase of R/bak® Cushion Mounting Materials and actual end use on the press. Periodically questions arise about proper storage of R/bak materials during this interval, specifically their shelf life and the best methods of storage to maintain product quality.



R/bak® SF, SS, and U – Corrugated Cushion Mounting Materials

Normal warehouse conditions generally are suitable for safe storage of R/bak products. R/bak rolls outside original containers should be stored away from prolonged exposure to UV or direct sunlight. Rotate inventory to avoid exceeding product shelf life.

When stored in these conditions, Rogers Corporation warrants a minimum shelf life of two years from the date of purchase for R/bak SF, SS and U corrugated products.

R/bak® SA 3000, SA 2000, SR 2000, PSA1, PSA2 and U1A Cushion Mounting Tapes

For protection against moisture, contamination, or handling damage, return unused R/bak tape to its original container for storage. Doing so will allow the roll to be easily identified.

When the tapes are outside of the original containers, avoid prolonged exposure to UV light and/or temperatures over 120°F (49°C).

When stored in these conditions, Rogers Corporation warrants a minimum shelf life of one year from the date of purchase for R/bak PSA1, PSA2, SR2000, SA2000, SA3000 and U1A mounting tapes.



The information contained in this Application Note is intended to assist you using Rogers R/bak Cushion Mounting Materials. It is not intended to and does not create any warranties, express or implied, including any warranty of merchantability or fitness for a particular purpose or that the results shown in this Application Note will be achieved by a user for a particular purpose. The user should determine the suitability of Rogers R/bak Cushion Mounting Materials for each application.

The world runs better with Rogers.®

The Rogers logo, The world runs better with Rogers., and R/bak are licensed trademarks of Rogers Corporation.
 © 2001, 2006, 2007 Rogers Corporation. All rights reserved. Printed in USA.
 7106-0807-PDF, Publication #12-046